

Product datasheet for **TP310536L**


OXLD1 (NM_001039842) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 17 open reading frame 90 (C17orf90), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210536 protein sequence Red =Cloning site Green =Tags(s) MLLRRVVEGGRAVAAAVRGSGARRFSSPDCCQRLPGGGSFLQRHHPGAQAPDGRRKFGTDHVEVGSQA GA DGTRPPKASLPPELQPPTNCCMSGCPNCVWEYADRLLQHFQDGGGERALAALEEHVADENLKAFLRMEI R LHTRCGG TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	15.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001034931</u>
Locus ID:	339229


[View online »](#)

UniProt ID: Q5BKU9
 RefSeq Size: 635
 Cytogenetics: 17q25.3
 RefSeq ORF: 441
 Synonyms: C17orf90

Product images:


Coomassie blue staining of purified OXLD1 protein (Cat# [TP310536]). The protein was produced from HEK293T cells transfected with OXLD1 cDNA clone (Cat# [RC210536]) using MegaTran 2.0 (Cat# [TT210002]).