

Product datasheet for **TP310562**

Acid Phosphatase 2 (ACP2) (NM_001610) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens acid phosphatase 2, lysosomal (ACP2), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210562 representing NM_001610 Red =Cloning site Green =Tags(s)
	<p>MAGKRXXGWSRAALLQLLLGVNLVMPPTRARSLRFVTLTYRHGDRSPVKTPKDPYQEEWPQGFGQL T KEGMLQHWELGQALRQRYHGFLNTSYHRQEVYVRSTDFDRTLMSAEANLAGLFPPNGMQRFNPNISW QPI PVHTVPITEDRLLKFPLGPCPRYEQLQNETRQTPEYQNESSRNAQFLDMVANETGLTDLTLETVWNVYDT LFCEQTHGLRLPPWASPQTMQRLSRLKDFSRFLFGIYQQAELARLQGGVLLAQIRKNLTLMATTSQLPK LLVYSAHDTTLVALQMALDVYNGEQAPYASCHIFELYQEDSGNFSVEMYFRNESDKAPWPLSLPGCPHRC PLQDFLRLTEPVVPKDWQQECQLASGPADTEVIVALAVCGSILFLLIVLLTLVFRMQAQPYPGRHVADG EDHA</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	45.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001601
Locus ID:	53
UniProt ID:	P11117
RefSeq Size:	2218
Cytogenetics:	11p11.2 11p12-p11
RefSeq ORF:	1272
Synonyms:	LAP
Summary:	The protein encoded by this gene belongs to the histidine acid phosphatase family, which hydrolyze orthophosphoric monoesters to alcohol and phosphate. This protein is localized to the lysosomal membrane, and is chemically and genetically distinct from the red cell acid phosphatase. Mice lacking this gene showed multiple defects, including bone structure alterations, lysosomal storage defects, and an increased tendency towards seizures. An enzymatically-inactive allele of this gene in mice showed severe growth retardation, hair-follicle abnormalities, and an ataxia-like phenotype. Alternatively spliced transcript variants have been found for this gene. A C-terminally extended isoform is also predicted to be produced by the use of an alternative in-frame translation termination codon via a stop codon readthrough mechanism. [provided by RefSeq, Oct 2017]
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Lysosome, Riboflavin metabolism

Product images:



Coomassie blue staining of purified ACP2 protein (Cat# TP310562). The protein was produced from HEK293T cells transfected with ACP2 cDNA clone (Cat# [RC210562]) using MegaTran 2.0 (Cat# [TT210002]).