

Product datasheet for **TP310581**

NUP62 (NM_016553) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human nucleoporin 62kDa (NUP62), transcript variant 2, 20 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC210581 protein sequence
Red=Cloning site **Green**=Tags(s)

MSGFNFGGTGAPTGGFTFGTAKTATTTTPATGFSFSTSGTGGFNFNGAPFQPATSTPSTGLFSLATQTPATQ
TTGFTFGTATLASGGTGFSGLGIGASKLNLNTAATPAMANPSGFLGSSNLTAISSTVTSSQGTAPTGF
VFGPSTTSVAPATTSGGFSFTGGSTAQPSGFNIGSAGNSAQPTAPATLPFTPATPAATTAGATQPAAPT
TATITSTGPSLFASIATAPTSSATTGLSLCTPVTTAGAPTAGTQGFSCLKAPGAASGTSTTTSTAATATAT
TTTSSSTTGFA LNLKPLAPAGIPSNTAAVTAPPGPAAAGAAASSAMTYAQLES LINKWSLELEDQERH
FLQQATQVNAWDRTLIENGEKITSLHREVEKVKLDQKRLDQELDFILSQQKELEDLLSPLEELVKEQSGT
IYLQHADEEREKTYKLAENIDAQLKRMAQDLKDIIHLNTSGAPADTSDPLQKICKILNAHMDSLQWIDQ
NSALLQRKVEEVTKVCEGRRKEQERSFRITFD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 53.1 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_057637](#)

Locus ID: 23636

UniProt ID: [P37198](#)

RefSeq Size: 3602

Cytogenetics: 19q13.33


RefSeq ORF: 1566

Synonyms: IBSN; p62; SNDI

Summary: The nuclear pore complex is a massive structure that extends across the nuclear envelope, forming a gateway that regulates the flow of macromolecules between the nucleus and the cytoplasm. Nucleoporins are the main components of the nuclear pore complex in eukaryotic cells. The protein encoded by this gene is a member of the FG-repeat containing nucleoporins and is localized to the nuclear pore central plug. This protein associates with the importin alpha/beta complex which is involved in the import of proteins containing nuclear localization signals. Multiple transcript variants of this gene encode a single protein isoform. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Transcription Factors

Product images:



Coomassie blue staining of purified NUP62 protein (Cat# TP310581). The protein was produced from HEK293T cells transfected with NUP62 cDNA clone (Cat# [RC210581]) using MegaTran 2.0 (Cat# [TT210002]).