

## Product datasheet for **TP310584**

### **MORN4 (NM\_178832) Human Recombinant Protein**


#### **Product data:**

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human MORN repeat containing 4 (MORN4), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210584 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MTLTKGSFTYSSGEEYRGWEKGRRHGFGQLMFADGGTYLGHFENGLFNGFGVLTFSDGSRYEGEFAQG K FNGVGVFIRYDNMTFEGEFKNGRVDGFGLLTFPDGSHGIPRNEGLFENNKLLRREKCSAIVQRAQSASKS ARNLTA  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	16.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_849154</a></u>
Locus ID:	118812


[View online »](#)

UniProt ID:	<u>Q8NDC4</u>
RefSeq Size:	2441
Cytogenetics:	10q24.2
RefSeq ORF:	438
Synonyms:	bA548K23.4; C10orf83; rtp; UTA
Summary:	Plays a role in promoting axonal degeneration following neuronal injury by toxic insult or trauma.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome

### Product images:



Coomassie blue staining of purified MORN4 protein (Cat# TP310584). The protein was produced from HEK293T cells transfected with MORN4 cDNA clone (Cat# [RC210584]) using MegaTran 2.0 (Cat# [TT210002]).