

Product datasheet for TP310589

MAT2B (NM_013283) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human methionine adenosyltransferase II, beta (MAT2B), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210589 protein sequence Red =Cloning site Green =Tags(s)

MVGREKELSIHFVPGSCRLVEEEVNIPNRRVLVTGATGLLGRAVHKEFQQNNWHAVGCGFRRARPKFEQV
NLLDSNAVHHIHDQPHVIVHCAAERRPDVVENQPDAASQLNVDASGNLAKEAAVGAFLIYSSDYVF
DGTNPPYREGDIPAPLNLYGKTKLDGEKAVLENNLGA AVLRIPILYGEVEKLEESAVTVMFDKVQFSNKS
ANMDHWQQRFPPTHVKDVATVCRQLAEKRMLDPSIKGTFHWSGNEQMTKYEMACAIADAFNLPSHLRPIT
DSPVLGAQRPRNAQLDCKLETLGIGQRTPFIRIGIKESLWPFLLIDKRWRQTVFH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	37.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_037415
Locus ID:	27430



[View online »](#)

UniProt ID:	Q9NZL9
RefSeq Size:	2154
Cytogenetics:	5q34
RefSeq ORF:	1002
Synonyms:	MAT-II; MATIIbeta; Nbla02999; SDR23E1; TGR
Summary:	The protein encoded by this gene belongs to the methionine adenosyltransferase (MAT) family. MAT catalyzes the biosynthesis of S-adenosylmethionine from methionine and ATP. This protein is the regulatory beta subunit of MAT. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Nov 2012]
Protein Pathways:	Cysteine and methionine metabolism, Metabolic pathways, Selenoamino acid metabolism

Product images:



Coomassie blue staining of purified MAT2B protein (Cat# TP310589). The protein was produced from HEK293T cells transfected with MAT2B cDNA clone (Cat# [RC210589]) using MegaTran 2.0 (Cat# [TT210002]).