

Product datasheet for **TP310598**

NDRG2 (NM_201541) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human NDRG family member 2 (NDRG2), transcript variant 8, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210598 protein sequence Red =Cloning site Green =Tags(s)

MAELQEVQITEEKPLLPGQTPEAAKTHSVETPYGSVTFVYGTGPKPKRPAILTYHDVGLNYKSCFQPLFQ
FEDMQEIIQNFVRVHVDAPGMEEGAPVFPLGYQYPSLDELADMIPCVLQYLNFTIIGVGVGAGAYILAR
YALNHPDTVEGLVLINIDPNAKGWMDWAAHKLTLTSSIPEMILGHLFSQEELSGNSELIQYRNIITHA
PNLDNIELYWNSYNNRRDLNFERGGDITLRCPVMLVVGDAQPHEDAVECNSKLDPTRTSFLKMADSGGQ
PQLTQPGKLTFAFKYFLQGMGYMASSCMTRLSRSRTASLTSAAASVDGNRSRRTLSQSSESGLSSGPPG
HTMEVSC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	39.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_963835</u>
Locus ID:	57447



[View online »](#)

UniProt ID: [Q9UN36](#)

RefSeq Size: 2077


Cytogenetics: 14q11.2

RefSeq ORF: 1071

Synonyms: SYLD

Summary: This gene is a member of the N-myc downregulated gene family which belongs to the alpha/beta hydrolase superfamily. The protein encoded by this gene is a cytoplasmic protein that may play a role in neurite outgrowth. This gene may be involved in glioblastoma carcinogenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2017]

Product images:



Coomassie blue staining of purified NDRG2 protein (Cat# TP310598). The protein was produced from HEK293T cells transfected with NDRG2 cDNA clone (Cat# [RC210598]) using MegaTran 2.0 (Cat# [TT210002]).