

Product datasheet for TP310628M

DDX4 (NM_024415) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human DEAD (Asp-Glu-Ala-Asp) box polypeptide 4 (DDX4), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210628 protein sequence Red =Cloning site Green =Tags(s)
	<p>MGDEDWEAEINPHMSSYVPIFEKDRYSGENGDNFNRTPASSEMDDGPSRRDHFMKSGFASGRNFGNRDA GECNKRDNTSTMGGFGVKGKSFNGRFGNSRFEDGDSSGFWRESSNDCEDNPTRNRGFSKRGGYRDGNNSE ASGPYRRGGRGSRGCRGGFGLGSPNNDLDPDECMQRTGGLFGSRRPVLSTGNGDTSQSRSGSGSERGG YKGLNEEVITGSGKNSWKSEAEGGESSDTQGPKVYIPIPPPEDEDSIFAHYQTGINFDKYDTILVEVSG HDAPPAILTFEEANLCQTLNNNIAKAGYTKLTPVQKYSIPIILAGRDLMACAQTGSGKTA AFLLPILAHM MHDGITASRFKELQEPECIIVAPTRELVNQIYLEARKFSFGTCVRAVVIYGGTQLGHSIRQIVQGCNILC ATPGR LMDIIGKEKIGLKQIKYLVLDEADRMLDMGFGPEMKKLISCPGMPSKEQRQTLMFSAFP EEIQR LAAEFLKSNYLFVAVGQVGGACRDVQQTVLQVQVQFSKREKLVEILRNIGDERTMV FVETKKKADFIATFL CQEKISTTSIHGDREQREREQALGDFRFGKCPVLVATSVAARGLDIENVQHVINF DLPSTIDEVYVHRIGR TGRCGNTGRAISFFDLESDNHLAQPLVKVLTDAQQDVP AWLEEIAFSTYIPGFSGSTRGNVFASVDTRKG KSTLNTAGFSSSQAPNPVDESWD</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	79.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_077726
Locus ID:	54514
UniProt ID:	Q9NQI0
RefSeq Size:	2880
Cytogenetics:	5q11.2
RefSeq ORF:	2172
Synonyms:	VASA
Summary:	DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene encodes a DEAD box protein, which is a homolog of VASA proteins in <i>Drosophila</i> and several other species. The gene is specifically expressed in the germ cell lineage in both sexes and functions in germ cell development. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009]

Product images:



Coomassie blue staining of purified DDX4 protein (Cat# [TP310628]). The protein was produced from HEK293T cells transfected with DDX4 cDNA clone (Cat# [RC210628]) using MegaTran 2.0 (Cat# [TT210002]).