

## Product datasheet for **TP310643**

### GFUS (NM\_003313) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tissue specific transplantation antigen P35B (TSTA3), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210643 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MGEPQGS MRILVTGG SGLVGKAIQKV VADGAGLPGEDWV FVSSKDADLTDTAQTRALFEKVQPTHVIHLA  
AMVGG LFRNIKYNLDFWRKNVHMNDNLHSAFEV GARKVVSCLSTCIFPDKTTYPIDETMIHNGPPHNSN  
FGYSYAKRMIDVQNRAYFQQYGCTFTAVIPTNVFGPHDNFNIEDGHVLPGLIHKVHLAKSSGSALT VWGT  
GNPRRQFIYSLDLAQLFIWVLR EYNEVEPIILSVGEEDEVSIKEAAEAVVEAMDFHGEVTFD TTKSDGQF  
KKTASNSKLR TYLPDFRFTP FKQAVKETCAWFTD NYEQARK

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	35.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_003304</a>
Locus ID:	7264



[View online »](#)

UniProt ID: [Q13630](#), [A0A140VKC8](#)

RefSeq Size: 1363

Cytogenetics: 8q24.3

RefSeq ORF: 963


Synonyms: FX; P35B; SDR4E1; TSTA3

**Summary:** Tissue specific transplantation antigen P35B is a NADP(H)-binding protein. It catalyze the two-step epimerase and the reductase reactions in GDP-D-mannose metabolism, converting GDP-4-keto-6-D-deoxymannose to GDP-L-fucose. GDP-L-fucose is the substrate of several fucosyltransferases involved in the expression of many glycoconjugates, including blood group ABH antigens and developmental adhesion antigens. Mutations in this gene may cause leukocyte adhesion deficiency, type II. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome

**Protein Pathways:** Amino sugar and nucleotide sugar metabolism, Fructose and mannose metabolism, Metabolic pathways

### Product images:



Coomassie blue staining of purified TSTA3 protein (Cat# TP310643). The protein was produced from HEK293T cells transfected with TSTA3 cDNA clone (Cat# [RC210643]) using MegaTran 2.0 (Cat# [TT210002]).