

Product datasheet for **TP310651M**

OIF (OGN) (NM_014057) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human osteoglycin (OGN), transcript variant 3, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210651 protein sequence Red =Cloning site Green =Tags(s)
	 MKTLLQSTLLLLLLVPLIKPAPPTQQDSRIIYDYGTDNFEEIFSQDYEDKYLDGKNIKEKETVIIPNEKS LQLQKDEAITPLPPKKENDEMPPTCLLCVCLSGSVYCEEVDIDAVPPLPKESAYLYARFNKIKKLTAKDFA DIPNLRRLDFTGNLIEDIEDGTFSKLSLLEELSLAENQLLKLPLVLPKLTFLNAKYNKIKSRGIKANAFK KLNNLTFLYLDHNALESVPLNLPESLRVIHLQFNNIASITDDTFCKANDTSYIRDRIIEIRLEGNPVILG KHPNSFICKRLPIGSYF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	31.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_054776</u>
Locus ID:	4969



[View online »](#)

UniProt ID: [P20774](#), [A8K0R3](#), [Q7Z532](#), [B4DI63](#)

RefSeq Size: 2786

Cytogenetics: 9q22.31


RefSeq ORF: 894

Synonyms: OG; OIF; SLRR3A

Summary: This gene encodes a member of the small leucine-rich proteoglycan (SLRP) family of proteins. The encoded protein induces ectopic bone formation in conjunction with transforming growth factor beta and may regulate osteoblast differentiation. High expression of the encoded protein may be associated with elevated heart left ventricular mass. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016]

Protein Families: Secreted Protein

Product images:



Coomassie blue staining of purified OGN protein (Cat# [TP310651]). The protein was produced from HEK293T cells transfected with OGN cDNA clone (Cat# [RC210651]) using MegaTran 2.0 (Cat# [TT210002]).