

Product datasheet for **TP310664M**

Serum Amyloid A (SAA1) (NM_000331) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human serum amyloid A1 (SAA1), transcript variant 1, 100 µg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC210664 protein sequence
Red=Cloning site **Green**=Tags(s)

MKLLTGLVFCSLVLGVSSRSFFSFLGEAFDGMWRAYSMDREANYIGSDKYFHARGNYDAAKRGPGGV
WAAEAISDARENIQRFFGHGAEDSLADQAANEWGRSGKDPNHFRPAGLPEKY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 13.4 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_000322](#)

Locus ID: 6288

UniProt ID: [P02735](#), [P0DJ18](#)

RefSeq Size: 678

Cytogenetics: 11p15.1




[View online »](#)

RefSeq ORF: 366

Synonyms: PIG4; SAA; SAA2; TP53I4

Summary: This gene encodes a member of the serum amyloid A family of apolipoproteins. The encoded preproprotein is proteolytically processed to generate the mature protein. This protein is a major acute phase protein that is highly expressed in response to inflammation and tissue injury. This protein also plays an important role in HDL metabolism and cholesterol homeostasis. High levels of this protein are associated with chronic inflammatory diseases including atherosclerosis, rheumatoid arthritis, Alzheimer's disease and Crohn's disease. This protein may also be a potential biomarker for certain tumors. Finally, antimicrobial activity against *S. aureus* and *E. coli* resides in the N-terminal portion of the mature protein. Alternate splicing results in multiple transcript variants that encode the same protein. A pseudogene of this gene is found on chromosome 11. [provided by RefSeq, Jul 2020]

Product images:



Coomassie blue staining of purified SAA1 protein (Cat# [TP310664]). The protein was produced from HEK293T cells transfected with SAA1 cDNA clone (Cat# [RC210664]) using MegaTran 2.0 (Cat# [TT210002]).