

Product datasheet for TP310670M

TSNAXIP1 (NM_018430) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human translin-associated factor X interacting protein 1 (TSNAXIP1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210670 protein sequence Red =Cloning site Green =Tags(s)

MGGHLSWPPTYTSGQTILQNRKPCSDDYRKRVGSCQQHPFRTAKPQYLEELENYLKELLLDLGTDSTQ
 ELRLQPYREIFFIEDFKTYKPLLSSIKNAYEGMLAHQREKIRALEPLKAKLVTNEDCNERILAMRAE
 EKYEISLLKKEKMNLLKLIDKKNEEKISLQSEVTKLRKNLAEEYLHYLSERDACKILIADLNELRYQRED
 MSLAQSPGIWGEDPVKLTALKMTRQDLTRTQMELNMMKANFGDVVPRRDFEMQEKTNKDLQEQLDT
 LRA
 SYEEVRKEHEILMQLHMSTLKERDQFFSELQEIQRTSTPRPDWTKCKDVGAGPERWQMLAEGKNSDQL
 V
 DVLLEEIGSGLLREKDFPGLGYGEAIPAFLRFDGLVENKKPSKKDWNLLKDAWKERLAEQKETFPDF
 FNFLEHRFGPSDAMAWAYTIFENIKIFHSNEVMSQFYAVLMGKRSENVYVTQKETVAQLLKEMTNADSQ
 NEGLLTMEQFNTVLKSTFPLKTEEQIQELMEAGGWHPPSSSNADLLNYRSLFMEDEEGQSEPFVQKLWEQ
 Y
 MDEKDEYLQQLKQELGIELHEEVTLPKLRGGLMTIDPSLDKQTVNTYMSQAFQLPESEMPEEGDEKEEAV
 VEILQTALERLQVIDIRRVGPREPEPAS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	76.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.


[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_060900
Locus ID:	55815
UniProt ID:	Q2TAA8
RefSeq Size:	2692
Cytogenetics:	16q22.1
RefSeq ORF:	1974
Synonyms:	TXI1
Summary:	Possible role in spermatogenesis.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified TSNAXIP1 protein (Cat# [TP310670]). The protein was produced from HEK293T cells transfected with TSNAXIP1 cDNA clone (Cat# [RC210670]) using MegaTran 2.0 (Cat# [TT210002]).