

Product datasheet for **TP310681L**

PRH2 (NM_005042) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human proline-rich protein HaellI subfamily 2 (PRH2), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210681 protein sequence Red =Cloning site Green =Tags(s) MLLILLSVALLAFSSAQDLDEDVSQEDVPLVISDGGDSEQFIDEERQGPPLGGQQSQPSAGDGNQDDGPQ QGPPQQGGQQQGPPPPQGKPGPPQQGGHPPPPQGRPQGPPQQGGHPRPPRGRPGPPQQGG HQQGPPP PPPGKPGPPPPQGGRPGPPQGQSPQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	16.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_005033</u>
Locus ID:	5555


[View online »](#)

UniProt ID:	<u>P02810</u>
RefSeq Size:	4323
Cytogenetics:	12p13.2
RefSeq ORF:	498
Synonyms:	db-s; pa; PIF-S; Pr; pr1/Pr2; PRH1; PRP-1/PRP-2
Summary:	This gene encodes a member of the heterogeneous family of proline-rich salivary glycoproteins. The encoded preproprotein undergoes proteolytic processing to generate one or more mature isoforms before secretion from the parotid and submandibular/sublingual glands. In western population this locus is commonly biallelic and encodes proline-rich protein (PRP) isoforms, PRP-1 and PRP-2. The reference genome encodes the PRP-1 allele. Certain alleles of this gene are associated with susceptibility to dental caries. This gene is located in a cluster of closely related salivary proline-rich proteins on chromosome 12. [provided by RefSeq, Oct 2015]
Protein Families:	Secreted Protein

Product images:



Coomassie blue staining of purified PRH2 protein (Cat# [TP310681]). The protein was produced from HEK293T cells transfected with PRH2 cDNA clone (Cat# [RC210681]) using MegaTran 2.0 (Cat# [TT210002]).