

Product datasheet for **TP310691**

H2A.Z (H2AFZ) (NM_002106) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant protein of human H2A histone family, member Z (H2AFZ), 20 µg |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >RC210691 protein sequence Red =Cloning site Green =Tags(s) |
| | MAGGKAGKDSGKAKTKAVSRSQRAGLQFPVGRIHRHLKSRRTTSHGRVGATAAVYSAAILEYLTAEVLELAGNASKDLKVKRITPRHLQLAIRGDEELDSLKATIAGGGVIPHHKSLIGKKGQKTV |
| | TRTRPLEQKLISEEDLAANDILDYKDDDDKV |
| Tag: | C-Myc/DDK |
| Predicted MW: | 13.4 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Preparation: | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps. |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_002097 |
| Locus ID: | 3015 |
| UniProt ID: | P0C0S5 |
| RefSeq Size: | 951 |
| Cytogenetics: | 4q23 |



[View online »](#)

RefSeq ORF: 384

Synonyms: H2A.z; H2A.Z-1; H2A/z; H2AFZ; H2AZ

Summary: Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene encodes a replication-independent member of the histone H2A family that is distinct from other members of the family. Studies in mice have shown that this particular histone is required for embryonic development and indicate that lack of functional histone H2A leads to embryonic lethality. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Systemic lupus erythematosus

Product images:



Coomassie blue staining of purified H2AFZ protein (Cat# TP310691). The protein was produced from HEK293T cells transfected with H2AFZ cDNA clone (Cat# [RC210691]) using MegaTran 2.0 (Cat# [TT210002]).