

## Product datasheet for **TP310762M**

### Apolipoprotein A I (APOA1) (NM\_000039) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human apolipoprotein A-I (APOA1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210762 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MKAAVLTLAVLFLTGSQARHFWQQDEPPQSPWDRVKDLATVVVDVLKDSGRDYVSQFEGSALGKQLNLKL  
LDNWDSVTSTFSKLREQLGPVTQEFWDNLEKETEGLRQEMSKDLEEVKAKVQPYLDDFQKKWQEEMELYR  
QKVEPLRAELQEGARQKLHELQEKLSPLGEEMRDRARAHVDALRTHLAPYSDELQRQLAARLEALKENGG  
ARLAEYHAKATEHLSTLSEKAKPALEDLRQGLLPVLESFKVSFLSALEEYTKKLNTQ

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	28.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_000030</a>
Locus ID:	335
UniProt ID:	<a href="#">P02647</a> , <a href="#">A0A024R3E3</a>



[View online »](#)

RefSeq Size: 897

Cytogenetics: 11q23.3

RefSeq ORF: 801

Synonyms: apo(a); HPALP2

**Summary:** This gene encodes apolipoprotein A-I, which is the major protein component of high density lipoprotein (HDL) in plasma. The encoded preproprotein is proteolytically processed to generate the mature protein, which promotes cholesterol efflux from tissues to the liver for excretion, and is a cofactor for lecithin cholesterolacyltransferase (LCAT), an enzyme responsible for the formation of most plasma cholesteryl esters. This gene is closely linked with two other apolipoprotein genes on chromosome 11. Defects in this gene are associated with HDL deficiencies, including Tangier disease, and with systemic non-neuropathic amyloidosis. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein. [provided by RefSeq, Dec 2015]

**Protein Families:** Druggable Genome, Secreted Protein

**Protein Pathways:** PPAR signaling pathway

### Product images:



Coomassie blue staining of purified APOA1 protein (Cat# [TP310762]). The protein was produced from HEK293T cells transfected with APOA1 cDNA clone (Cat# [RC210762]) using MegaTran 2.0 (Cat# [TT210002]).