

Product datasheet for **TP310767**

DCK (NM_000788) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human deoxycytidine kinase (DCK), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210767 protein sequence Red =Cloning site Green =Tags(s)
	 MATPPKRSCPSFSASSEGTRIKKISIEGNI AAGKSTFVNILKQLCEDWEVWPEPVARWCNVQSTQDEFEE LTMSQKNGGNVLQMMYEKPERWSFTFQTYACLSRIRAQLASLNGKLDKDAEKPVLFERSVYSDRYIFASN LYESEC MNETEWTIYQDWHDMNMQFGQSLELDGIIYLQATPETCLHRIYLRGRNEEQGIPLEYLEKLHY KHESWLLHRTLKTNFDY LQEVPI LLDVNE DFKDKYESLVEKVKEFLSTL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	30.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000779
Locus ID:	1633
UniProt ID:	P27707



[View online »](#)

RefSeq Size: 2618

Cytogenetics: 4q13.3


RefSeq ORF: 780

Summary: Deoxycytidine kinase (DCK) is required for the phosphorylation of several deoxyribonucleosides and their nucleoside analogs. Deficiency of DCK is associated with resistance to antiviral and anticancer chemotherapeutic agents. Conversely, increased deoxycytidine kinase activity is associated with increased activation of these compounds to cytotoxic nucleoside triphosphate derivatives. DCK is clinically important because of its relationship to drug resistance and sensitivity. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Purine metabolism, Pyrimidine metabolism

Product images:



Coomassie blue staining of purified DCK protein (Cat# TP310767). The protein was produced from HEK293T cells transfected with DCK cDNA clone (Cat# [RC210767]) using MegaTran 2.0 (Cat# [TT210002]).