

Product datasheet for **TP310785**

C14orf177 (NM_182560) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 14 open reading frame 177 (C14orf177), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210785 protein sequence Red =Cloning site Green =Tags(s)
	MHRKEPGARLEATRGAARPHKQGTPMITRPSVSQLGEGKCPSSQHLQSLRHNKQHALTLTKARCCGECS TCFCTEEKSECQRHEETSPGSCNHQIMSASTISAFCATPRFKQLFKGTVEQMSQM
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	13.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_872366
Locus ID:	283598
RefSeq Size:	1289
Cytogenetics:	14q32.2
RefSeq ORF:	375



[View online »](#)

Synonyms: FLJ25773

Product images:



Coomassie blue staining of purified C14orf177 protein (Cat# TP310785). The protein was produced from HEK293T cells transfected with C14orf177 cDNA clone (Cat# [RC210785]) using MegaTran 2.0 (Cat# [TT210002]).