

Product datasheet for **TP310831**

MST1 (STK4) (NM_006282) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human serine/threonine kinase 4 (STK4), 20 µg

Species: Human

Expression Host: HEK293T

Expression cDNA >RC210831 protein sequence

Clone or AA Sequence: **Red**=Cloning site **Green**=Tags(s)

METVQLRNPPRRQLKKLDEDSLTKQPEEVFDVLEKLGEGSVYKAIHKETGQIVAIAIKQVPVESDLQEI
IKEISIMQQCDSPHWKYYGSYFKNTDLWIVMEYCGAGSVSDIIRLRNKTLTDEIATILQSTLKGLEYL
HFMRKIHRDIKAGNILLNTEGHAKLADFGVAGQLTDTMAKRNTVIGTPFWMAPEVIQEIGYNCVADIWSL
GITAIEMAEGKPPYADIHPMRAIFMIPTNPPPTFRKPELWSDNFTDFVKQCLVKSPEQRATATQLLQHPF
VRSKAGVSILRDLINEAMDVKLKRQESQQREVDQDDEENSEEDEMDSGTMVRAVGDDEMGTVRVASTMTDG
ANTMIEHDDTLPSQLGTMVINAEDDEEEEGTMKRRDET MQPAKPSFLEYFEQKEKENQINSFGKSVPGPLK
NSSDWKIPQDGDYEFLLKSWTVEDLQKRLALDPMMEQEIEEIRQKYQSKRQPILDAIEAKRRRQQNF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 55.4 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


RefSeq: [NP_006273](#)



[View online »](#)

Locus ID:	6789
UniProt ID:	Q13043
RefSeq Size:	6344
Cytogenetics:	20q13.12
RefSeq ORF:	1461
Synonyms:	KRS2; MST1; YSK3
Summary:	The protein encoded by this gene is a cytoplasmic kinase that is structurally similar to the yeast Ste20p kinase, which acts upstream of the stress-induced mitogen-activated protein kinase cascade. The encoded protein can phosphorylate myelin basic protein and undergoes autophosphorylation. A caspase-cleaved fragment of the encoded protein has been shown to be capable of phosphorylating histone H2B. The particular phosphorylation catalyzed by this protein has been correlated with apoptosis, and it's possible that this protein induces the chromatin condensation observed in this process. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	MAPK signaling pathway, Non-small cell lung cancer, Pathways in cancer

Product images:



Coomassie blue staining of purified STK4 protein (Cat# TP310831). The protein was produced from HEK293T cells transfected with STK4 cDNA clone (Cat# [RC210831]) using MegaTran 2.0 (Cat# [TT210002]).