

Product datasheet for **TP310842L**


CCDC48 (EFCC1) (NM_024768) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 48 (CCDC48), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210842 representing NM_024768 Red =Cloning site Green =Tags(s) MTYFGHFGGANHAHTLGELEACIAMLVEQLRTQGCGGRTLGTSEEEAELQQKVEENEHLRLELQMVETER VRSLLEEKLVDVLQLLQRLRDLNISKRALGKILLSTLDAFRDPTHEGRPSAAILDALHQAACQLLR RQPSAPASAAAALTNPLLVS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	65.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_079044
Locus ID:	79825
UniProt ID:	Q9HA90
RefSeq Size:	1963


[View online »](#)

Cytogenetics: 3q21.3
 RefSeq ORF: 483
 Synonyms: C3orf73; CCDC48

Product images:


Coomassie blue staining of purified EFCC1 protein (Cat# [TP310842]). The protein was produced from HEK293T cells transfected with EFCC1 cDNA clone (Cat# [RC210842]) using MegaTran 2.0 (Cat# [TT210002]).