

Product datasheet for **TP310870M**

CD200R (CD200R1) (NM_170780) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human CD200 receptor 1 (CD200R1), transcript variant 4, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210870 protein sequence Red =Cloning site Green =Tags(s)
	 MLCPWRTANLGLLLILTIFLVAASSSLCMDEKQITQNYSKVLAEVNTSWPVKMATNAVLCCPPIALRNLI IITWEIILRGQPSTKAYKKETNETKTNCTDERITWVSRPDQNSDLQIRTVAITHDGYRCIMVTPDGN FHRGYHLQVLVTPVTLFQNRNRTAVCKAVAGKPAAHISWIPEGDCATKQEYWSNGTVTVKSTCHWEVHN VSTVTCHVSHLTGNKSLYIELLPVPGAKKSAKLYIPYIILTIILTVGFIWLLKVNCGCRKYKLNKTEST PVVEEDEMOPYASYTEKNNPLYDTTNKVKASEALQSEVDTDLHTL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	36.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_740750</u>
Locus ID:	131450




[View online »](#)

UniProt ID: [Q8TD46](#)
RefSeq Size: 2203
Cytogenetics: 3q13.2
RefSeq ORF: 975
Synonyms: CD200R; HCRTR2; MOX2R; OX2R

Summary: This gene encodes a receptor for the OX-2 membrane glycoprotein. Both the receptor and substrate are cell surface glycoproteins containing two immunoglobulin-like domains. This receptor is restricted to the surfaces of myeloid lineage cells and the receptor-substrate interaction may function as a myeloid downregulatory signal. Mouse studies of a related gene suggest that this interaction may control myeloid function in a tissue-specific manner. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Transmembrane

Product images:



Coomassie blue staining of purified CD200R1 protein (Cat# [TP310870]). The protein was produced from HEK293T cells transfected with CD200R1 cDNA clone (Cat# [RC210870]) using MegaTran 2.0 (Cat# [TT210002]).