

Product datasheet for **TP310871M**

LCN15 (NM_203347) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human lipocalin 15 (LCN15), 100 µg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC210871 protein sequence
Red=Cloning site **Green**=Tags(s)

MMSFLLGAILTLLWAPTAQAEVLLQPDFNAEKFSGLWYVSMASDCRVFLGKKDHLMSSTRAIRPTEEGG
LHVHMEFPGADGCNQVDAEYLKVGSEGHFRVPALGYLDVRIVDTDYSSFAVLIYKELEGALSTMVQLYS
RTQDVSPQALKSFQDFYPTLGLPKDMMVMLPQSDACNPESKEAP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 20.3 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_976222](#)

Locus ID: 389812


UniProt ID: [Q6UWW0](#)

RefSeq Size: 762



[View online »](#)

Cytogenetics: 9q34.3
RefSeq ORF: 552
Synonyms: PRO6093; UNQ2541

Product images:

Coomassie blue staining of purified LCN15 protein (Cat# [TP310871]). The protein was produced from HEK293T cells transfected with LCN15 cDNA clone (Cat# [RC210871]) using MegaTran 2.0 (Cat# [TT210002]).