

Product datasheet for **TP310873**


LOC283392 (NM_001025457) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hypothetical protein LOC283392 (LOC283392), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210873 representing NM_001025457 Red =Cloning site Green =Tags(s)
	MLPPPLRACGSRLSLPCAALSPRFPREPSNRFSAREKDGSGGKSEGVEREVRAFAFSPHSSGGVRRSSGWN VLA AA VTAATAAREAPEPYTLQLERKGWPRWERLE
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	11.8
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001020628
Locus ID:	283392
RefSeq Size:	430
Cytogenetics:	12q21.1
RefSeq ORF:	312



[View online »](#)

Product images:

Coomassie blue staining of purified LOC283392 protein (Cat# TP310873). The protein was produced from HEK293T cells transfected with LOC283392 cDNA clone (Cat# [RC210873]) using MegaTran 2.0 (Cat# [TT210002]).