

Product datasheet for **TP310877**

RAB11FIP4 (NM_032932) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RAB11 family interacting protein 4 (class II) (RAB11FIP4), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210877 protein sequence Red =Cloning site Green =Tags(s)

MAGGAGWSGAPAALLRSVRRRLREVFVCGRDPDGFLRVERVAALGLRFGQGEEVEKLVKYLDPNDLGRIN
FKDFCRGVFAMKGCCELLKDVLVSAGTLPCAPEIPDCVEQGSEVTGPTFADGELIPREPGFFPEDEEE
AMTLAPPEGPQELYTDSPMESTQSLEGSVSPAEDGGLGGLFLPEDKSLVHTPSMTTSDLSTHSTSLI
SNEEQFEDYGEGDDVDCAPSSPCPDDTRTNVYSDLGSSVSSAGQTPRKMHRVYNSSELLDVYCSQCCKK
INLLNDLEARLKNLKANSPNRKISSTAFGRQLMHSSNFSSNGSTEDLFRDSIDSCDNDITEKVSFLEKK
VTELENDSLTNGDLKSKLKQENTQLVHRVHELEEMVKDQETTAEQALEEEARRHREAYGKLEREKATEVE
LLNARVQQLLEEENELRTTVTRLKSQTEKLDEERQRMSDRLEDTSRLKDEMDLYKRMMDKLRQNRLEFQ
KEREATQELIEDLRKELEHLQMYKLD CERPGRGRSASSGLGEFNARAREVELEHEVKRLKQENYKLRDQN
DDLNGQILSLSLYEAKNFLAAQTKAQS LA AEIDTASRDELMEALKEQEEINFRLRQYMDKIILAILDHNP
SILEIKH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	71.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_116321](#)

Locus ID: 84440

UniProt ID: [Q86YS3](#)

RefSeq Size: 8665

Cytogenetics: 17q11.2


RefSeq ORF: 1911

Synonyms: FIP4-Rab11; RAB11-FIP4

Summary: The protein encoded by this gene interacts with RAB11 and is thought to be involved in bringing recycling endosome membranes to the cleavage furrow in late cytokinesis. Hypoxic conditions can lead to an upregulation of the encoded protein and enhance the metastatic potential of hepatocellular carcinoma. [provided by RefSeq, Oct 2016]

Protein Pathways: Endocytosis

Product images:



Coomassie blue staining of purified RAB11FIP4 protein (Cat# TP310877). The protein was produced from HEK293T cells transfected with RAB11FIP4 cDNA clone (Cat# [RC210877]) using MegaTran 2.0 (Cat# [TT210002]).