

Product datasheet for **TP310899L**

Epsin 2 (EPN2) (NM_148921) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human epsin 2 (EPN2), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210899 protein sequence Red =Cloning site Green =Tags(s)
	<p>MTTSSIRRQMKNI VNNYSEAEIKVREATSNDPWGPSSSLMTEIADLTYNVVFSEIMSMVWKRLNDHGKN WRHVYKAL TLLDYLIK TGSERVAQQCRENIFAIQTLKDFQYIDRDGKDQGINVREKSKQLVALLKDEERL KAERAQALKTKERMAQVATGMGSNQITFGRGSSQPNLSTSHSEQEYKGAGGSPASYHGSTSPRVSSELEQ ARPQTSGEELQLQLALAMSREVAEQEERLRRGDDLRLQMALEESRRD TVKIPKKKEHGSLPQQTTLLDL MDALPSSGPAQAEPWGPSASTNQTNPWGGAAPASTSDPWPSFGTKPAASIDPWGVPTGATAQSVPKN SDPWAASQQPASSAGKRASDAWGAVSTTKPVSVSGSFELFSNLNGTIKDDFSEFDNLRTSKKTAESVTSL PSQNNGTTSPDPFESQPLTVASSKPSSARKTPESFLGPNAALVNLDSLVRPAPPAQSLNPFLAPGAPAT SAPVNPVFQVNPQPLTLNQLRGSPVLGTSTSFPGPGVESMAVASMTS AAPQPALGATGSSLTPLGPAMM NMVGSVGI PPSAAQATGTTNPFL</p> <p>SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	62.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_683723
Locus ID:	22905
UniProt ID:	O95208
RefSeq Size:	4664
Cytogenetics:	17p11.2
RefSeq ORF:	1752
Synonyms:	EHB21
Summary:	This gene encodes a protein which interacts with clathrin and adaptor-related protein complex 2, alpha 1 subunit. The protein is found in a brain-derived clathrin-coated vesicle fraction and localizes to the peri-Golgi region and the cell periphery. The protein is thought to be involved in clathrin-mediated endocytosis. Alternate splicing of this gene results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]
Protein Pathways:	Endocytosis

Product images:



Coomassie blue staining of purified EPN2 protein (Cat# [TP310899]). The protein was produced from HEK293T cells transfected with EPN2 cDNA clone (Cat# [RC210899]) using MegaTran 2.0 (Cat# [TT210002]).