

Product datasheet for **TP310930M**

TOR1AIP2 (NM_145034) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human torsin A interacting protein 2 (TOR1AIP2), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210930 protein sequence Red =Cloning site Green =Tags(s)

MADSGLEPQEDSQKDLNDPSVNSQAQETIIASNAEEAEILHSACGLSKDHQEVETEGPESADTGDKS
ESPDEANVGKHPKDKTEDENKQSFLDGGKGGHLPSENLGKEPLDPDPSPDKVGRADAHLGSSVALP
KEASDGTGASQEPPTTDSQEAQSPGHSSAGQEGEDTLRRRLLAPEAGSHPQQTQKLEEIKENAQDTMRQI
NKKGFWSYGPVILVVLVAVVASSVNSYSSPAQQVPKNPALEAFLAQFSQLEDKFPQGSSFLWQRGRKF
LQKHLNASNPTEPATIIFTAAREGRETCLKLSHHVADAYTSSQKVSPIQIDGAGRTWQSDTVKLLVDLE
LSYGFENGQKAAVHHFESFPAGSTLIFYKYCDHENAAFKDVALVLTVLLEETLEASVGPREETEEKVRD
LLWAKFTNSDTPTSFNHMDSDKLSGLWSRISHLVLPVQPVSIEEQGLF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	51.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_659471</u>



[View online »](#)

Locus ID: 163590

UniProt ID: [Q8NFQ8](#)

RefSeq Size: 7912

Cytogenetics: 1q25.2


RefSeq ORF: 1410

Synonyms: IFRG15; LULL1; NET9

Summary: One of the two protein isoforms encoded by this gene is a type II integral membrane protein found in the endoplasmic reticulum (ER). The encoded protein is a cofactor for the ATPase TorsinA, regulating the amount of TorsinA present in the ER compared to that found in the nuclear envelope. Defects in this protein are a cause of early onset primary dystonia, a neuromuscular disease. The other isoform encoded by this gene is an interferon alpha responsive protein whose cellular role has yet to be determined. [provided by RefSeq, Mar 2017]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified TOR1AIP2 protein (Cat# [TP310930]). The protein was produced from HEK293T cells transfected with TOR1AIP2 cDNA clone (Cat# [RC210930]) using MegaTran 2.0 (Cat# [TT210002]).