

Product datasheet for TP310936M

GRASP1 (GRIPAP1) (NM_020137) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human GRIP1 associated protein 1 (GRIPAP1), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210936 protein sequence Red =Cloning site Green =Tags(s)

MAQALSEEEFQRMQAQLLELRTNNYQLSDELKNGVELTSLRQKVAYLDKEFSKAQKALSKSKKAQEVEV
LLENEMLQAKLHSQEEDFRLQNSTLMAEFSKLCSQMEQLEQENQQLKEGAAGAGVAQAGPLVDGELLRL
QAENTALQKNVAALQERYGKEAGKFSVSEGQDPPGGPAPTVLAPMPLAEVELKWEMEKEEKRLWEQL
QGLESSKQAETSRLQEELAKLSEKLEKQESFCRLQTEKETLFNDSRNKIEELQQRKEADHKAQLARTQK
LQQELEAANQSLAELRDQRQGERLEHAAALRALQDQVSIQSADAQEQQVEGLLAENNALRTSLAALEQIQT
AKTQELNMLREQTTGLAAELQQQAEYEDLMGQKDDLNSQLQESLRANSRLLQLEIGQEKEQLTQELQ
EARKSAEKRKAMLDELAMETLQEKSHKEELGAVRLRHEKEVLGVRARYERELRELHEDKKRQEEELRGQ
IREEKARTRELETLQQTVEELQAQVHSMGAKGWFERRLKEAEESSLQQQQQEQQEALKQCREQHAAELKG
KEEELQDVRDQLEQAQEERDCHLKTISLQEVKDTVDGQRILEKKGSAALKDLKRQLHLERKRADKLQE
RLQDILTNSKSRSGLEELVLESEMNSPSTQTGDSSSISSFSYREILREKESSAVPARSLSSSPQAQPPRP
AELSDEEVAELFQRLAETQEKWMLEEKVKHLEVSSASMAEDLCRKSIIETVYVMDSRIDVVAAGHTDR
SGLGSVLRDLVKPGDENLREMNNKLLQNMLEEQLTKNMHLHKDMEVLSQEIVRLSKECVGPPDPDLEPGET
S

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	95.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_064522](#)

Locus ID: 56850

UniProt ID: [Q4V328](#)

RefSeq Size: 3052

Cytogenetics: Xp11.23


RefSeq ORF: 2523

Synonyms: GRASP-1

Summary: This gene encodes a guanine nucleotide exchange factor for the Ras family of small G proteins (RasGEF). The encoded protein interacts in a complex with glutamate receptor interacting protein 1 (GRIP1) and plays a role in the regulation of AMPA receptor function. [provided by RefSeq, Aug 2013]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified GRIPAP1 protein (Cat# [TP310936]). The protein was produced from HEK293T cells transfected with GRIPAP1 cDNA clone (Cat# [RC210936]) using MegaTran 2.0 (Cat# [TT210002]).