

Product datasheet for **TP310940L**

TPRG1 (NM_198485) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tumor protein p63 regulated 1 (TPRG1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210940 protein sequence Red =Cloning site Green =Tags(s)
	<p>MSTIGSFEGFQAVSLKQEGDDQPSETDHLSEEDPMPRQISRQSSVTESTLYPNPYHQPYISRKYFATR PGAJETAMEDLKGHVAETSGETIQFWLLTKIDHWNNEKERILLVTDKTLICKYDFIMLSCVQLQRIPL SAVYRICLGKFTFPGMSLDKRQGEGLRIYWGSPPEEQSLLSRWNPWSTVPYATFTEHPMKYTSEKFLEIC KLSGFMSKLVPAIQNAHKNSTGSGRGKKLMVLTEPILIETYTGLMFIGNRNKLGYSLARGSIGF</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	31 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_940887</u>
Locus ID:	285386
UniProt ID:	<u>Q6ZUI0</u>



[View online »](#)

RefSeq Size: 3665

Cytogenetics: 3q28

RefSeq ORF: 825

Synonyms: FAM79B

Product images:


Coomassie blue staining of purified TPRG1 protein (Cat# [TP310940]). The protein was produced from HEK293T cells transfected with TPRG1 cDNA clone (Cat# [RC210940]) using MegaTran 2.0 (Cat# [TT210002]).