

Product datasheet for TP310999

LOC285033 (NM_001037228) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hypothetical protein LOC285033 (LOC285033), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC210999 protein sequence Red =Cloning site Green =Tags(s) MKSVYNIRPLTPQYVRLFLVAGNFFRIQKLLMEASVFIKEGAVCGIARPKTSRVNSSQDQIQVASENTH SGLHQRPASGARLPASCGQELPLFSQRAVGGLPRKQSEAHTPASFLPCS TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	13 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001032305</u>
Locus ID:	285033
RefSeq Size:	537
Cytogenetics:	2q11.2
RefSeq ORF:	363


[View online »](#)

Synonyms: MGC126779; MGC126781

Product images:



Coomassie blue staining of purified LOC285033 protein (Cat# TP310999). The protein was produced from HEK293T cells transfected with LOC285033 cDNA clone (Cat# [RC210999]) using MegaTran 2.0 (Cat# [TT210002]).