

## Product datasheet for TP311013L

### B Raf (BRAF) (NM\_004333) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human v-raf murine sarcoma viral oncogene homolog B1 (BRAF), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211013 protein sequence Red=Cloning site Green=Tags(s)

MAALSGGGGGGAEPGQALFNGDMEPEAGAGAGAAAASSAADPAIPEEVWNIKQMIKLTQEHEALLDKFGG  
EHNPPSIYLEAYEYTSKLDALQQREQLLESLGNGTDFSVSSSASMDTVTSSSSSSSLVLPSSLSVFQN  
PTDVARSNPKSPQKPIVRVFLPNKQRTVVPARCGVTVRDSLKKALMMRGLIPECCAVYRIQDGEKKPIGW  
DTDISWLTGEELHVEVLENVPLTTHNFVRKTFFTLAFCDFCRKLLFQGFRCQTCGYKFHQRCEVPLMC  
VNYDQLDLLFVSKFFEHHPIPQEEASLAETALTSGSSSPASDSIGPQILTSPSPSKSIPIQPFRPAD  
EDHRNQFGQRDRSSAPNVHINTIEPVNIDDLIRDQGFRGDGGSTTGLSATPPASLPGSLTNVKALQKSP  
GPQRERKSSSSSEDRNRMKTLGRRDSSDDWEIPDQGQITVGQRIGSGSFGTVYKKGWHGDVAVKMLNVTAP  
TPQQQLQAFKNEVGVLKTRHVNILLFMGYSTKPQLAIVTQWCEGSSLYHHLHIETKFEMIKLIDIARQT  
AQGMDYLHAKSIIHRDLKSNINFLHEDLTVKIGDFGLATVKSRSWGSQHQFEQLSGSILWMAPEVIRMQDK  
NPYSFQSDVYAFGIVLYELMTGQLPYSNINNRDQIIFMVGGRGYLSPDLSKVRSNCPKAMKRLMAECLKKK  
RDERPLFPQILASIELLARSPLKIHRSASEPSLNRAFQTEDFSLYACASPKTPIQAGGYGAFPVH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



Tag:	C-Myc/DDK
Predicted MW:	84.3 kDa
Concentration:	>0.1 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	BRAF kinase activity was measured in an HTRF® assay. Varying concentrations of BRAF were added to a reaction mix containing ATP and a biotinylated kinase substrate (HTRF substrate 2) and was incubated at 37C for phosphorylation. HTRF detection reagents were then added, the reaction was incubated for 30 minutes at room temperature. Time-resolved fluorescent signal (Delta R) was measured on a Flexstation 3 microplate reader.



[View online »](#)

<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_004324</a>
<b>Locus ID:</b>	673
<b>UniProt ID:</b>	<a href="#">P15056</a>
<b>RefSeq Size:</b>	2949
<b>Cytogenetics:</b>	7q34
<b>RefSeq ORF:</b>	2298
<b>Synonyms:</b>	B-raf; B-RAF1; BRAF1; NS7; RAFB1
<b>Summary:</b>	<p>This gene encodes a protein belonging to the RAF family of serine/threonine protein kinases. This protein plays a role in regulating the MAP kinase/ERK signaling pathway, which affects cell division, differentiation, and secretion. Mutations in this gene, most commonly the V600E mutation, are the most frequently identified cancer-causing mutations in melanoma, and have been identified in various other cancers as well, including non-Hodgkin lymphoma, colorectal cancer, thyroid carcinoma, non-small cell lung carcinoma, hairy cell leukemia and adenocarcinoma of lung. Mutations in this gene are also associated with cardiofaciocutaneous, Noonan, and Costello syndromes, which exhibit overlapping phenotypes. A pseudogene of this gene has been identified on the X chromosome. [provided by RefSeq, Aug 2017]</p>
<b>Protein Families:</b>	Druggable Genome, Protein Kinase
<b>Protein Pathways:</b>	Acute myeloid leukemia, Bladder cancer, Chemokine signaling pathway, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Focal adhesion, Glioma, Insulin signaling pathway, Long-term depression, Long-term potentiation, MAPK signaling pathway, Melanoma, mTOR signaling pathway, Natural killer cell mediated cytotoxicity, Neurotrophin signaling pathway, Non-small cell lung cancer, Pancreatic cancer, Pathways in cancer, Progesterone-mediated oocyte maturation, Prostate cancer, Regulation of actin cytoskeleton, Renal cell carcinoma, Thyroid cancer, Vascular smooth muscle contraction

Product images:



Coomassie blue staining of purified BRAF protein (Cat# [TP311013]). The protein was produced from HEK293T cells transfected with BRAF cDNA clone (Cat# [RC211013]) using MegaTran 2.0 (Cat# [TT210002]).