

Product datasheet for **TP311051L**

C2orf68 (NM_001013649) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 2 open reading frame 68 (C2orf68), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211051 protein sequence Red =Cloning site Green =Tags(s) MEAGPHPRPGHCCKPGGR LDMNHGFVHHIRRNQIARDDYDKVKQA AKEKVR RRHTPAPTRPRKPD L QVY LPRHRDVSAHPRNP DYEESSSGGSELEPSGHQLFCLEYEADSGEVT SVIVYQGDDPGKVSEK VSAH TPLDPPMREALKLRIQEEIAKRQSQH TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	18.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001013671</u>
Locus ID:	388969
UniProt ID:	<u>Q2NWX9</u>



[View online »](#)

RefSeq Size: 4283

Cytogenetics: 2p11.2

RefSeq ORF: 498

Synonyms: HCRCN81

Product images:


Coomassie blue staining of purified C2orf68 protein (Cat# [TP311051]). The protein was produced from HEK293T cells transfected with C2orf68 cDNA clone (Cat# [RC211051]) using MegaTran 2.0 (Cat# [TT210002]).