

## Product datasheet for **TP311055L**

### **C6ORF154 (LRRC73) (NM\_001012974) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 6 open reading frame 154 (C6orf154), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211055 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MLPSSIQISGEPLSGAEVRDICRGLRDNAVRLLSLRGCRLCDRDFGRICRALAGATSLAQLNINLGVSS  
PSRIKQLAEALRTNRSIQSLFLHGSPPLTDAGLALLNPALALHPALVALDLGDCMLGDEAINLICGLLPPD  
GAKSGLKELTLANPGITPKGWSRLAIAVAHSSQVRVNLNDYNPLGDHVAGMLAVAVASSRTLEVLDLEG  
TGLTNQSAQTLDMVENYPTALRSLVLAENSISPELQQQICDLLSEGEVEEVAGGAGDTQEWERGREPA  
AHQRGSSSWMCPSPSSQMVLMTSGLGDSLLAETEM


**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	33.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001012992</a></u>
Locus ID:	221424



[View online »](#)

UniProt ID: [Q5JTD7](#)  
RefSeq Size: 1805  
Cytogenetics: 6p21.1  
RefSeq ORF: 948  
Synonyms: C6orf154

**Product images:**

Coomassie blue staining of purified LRRC73 protein (Cat# [TP311055]). The protein was produced from HEK293T cells transfected with LRRC73 cDNA clone (Cat# [RC211055]) using MegaTran 2.0 (Cat# [TT210002]).