

Product datasheet for TP311067L

CYorf15A (NM_001005852) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome Y open reading frame 15A (CYorf15A), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211067 protein sequence Red =Cloning site Green =Tags(s) MEEAGLCGLREKADMLCNSSEHDILQHQDSNCSATSNKHLLDEEGRDFITKNRSWVSPVHCTQESRREL PEQEVAPPSGQQALQCNRNKEKVLERQGLSLSPRLERGGTFMGHGSLKLPVLVLSLQTLT TR TRPLE QKLISEEDLA AND ILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	14.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001005852
Locus ID:	246126
RefSeq Size:	872
Cytogenetics:	Yq11.222-q11.223
RefSeq ORF:	393


[View online »](#)

Synonyms: FLJ33216; MGC131732

Product images:



Coomassie blue staining of purified CYorf15A protein (Cat# [TP311067]). The protein was produced from HEK293T cells transfected with CYorf15A cDNA clone (Cat# [RC211067]) using MegaTran 2.0 (Cat# [TT210002]).