

## Product datasheet for **TP311076M**

### UBL5 (NM\_001048241) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ubiquitin-like 5 (UBL5), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211076 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MIEVVCNDRLGKKVRVKCNTDDTIGDLKLLIAAQTGTRWNKIVLKKWYTIKFDHVS LGDYEIHDGMNLEL YYQ
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	8.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001041706</a>
Locus ID:	59286
UniProt ID:	<a href="#">Q9BZL1</a> , <a href="#">A0A024R7B0</a>
RefSeq Size:	469
Cytogenetics:	19p13.2



[View online »](#)

RefSeq ORF: 219

Synonyms: HUB1

**Summary:** This gene encodes a member of a group of proteins similar to ubiquitin. The encoded protein is not thought to degrade proteins like ubiquitin but to affect their function through being bound to target proteins by an isopeptide bond. The gene product has been studied as a link to predisposition to obesity based on its expression in *Psammomys obesus*, the fat sand rat, which is an animal model for obesity studies. Variation in this gene was found to be significantly associated with some metabolic traits (PMID: 15331561) but not associated with childhood obesity (PMID: 19189687). Pseudogenes of this gene are located on chromosomes 3, 5 and 17. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jan 2013]

### Product images:



Coomassie blue staining of purified UBL5 protein (Cat# [TP311076]). The protein was produced from HEK293T cells transfected with UBL5 cDNA clone (Cat# [RC211076]) using MegaTran 2.0 (Cat# [TT210002]).