

Product datasheet for TP311163L

PTPRS (NM_130853) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protein tyrosine phosphatase, receptor type, S (PTPRS), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC211163 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)


MAPTWGPGMVSVWGPMLLVLLVGGCAAEEPPRFIPEKPKDQIGVSGGVASFVCQATGDPKPRVTWNKKG
KKVNSQRFETIEFDESAGAVLRIQPLRTPRDENVYECVAQNSVGEITVHAKLTVLREDQLPSGFPNIDMG
PQLKVERTRTATMLCAASGNPDPEITWFKDFLPVDPSASNGRIKQLRSGALQIESSEETDQGGYECVAT
NSAGVRYSSPANLYVRVRRVAPRFSILPMSHEIMPGGNVNITCVAVGSPMPYVKWMQGAEDLTPEDDMPV
GRNVLELTDVKDSANYTCVAMSSLGVIEAVAQITVKSLPKAPGTPMVTENTATSITITWDSGNPDVSY
VIEYKSKSQDGPYQIKEDITTRYSIGGLSPNSEYEIWWSAVNSIGQGPPSESVVTRTGEQAPASAPRNV
QARMLSATTMIVQWEEPVEPNLIRGYRVYITMEPEHPVGNWQKHNVDDSLTTVGSLLLEDETYTVRVL
FTSVGDGPLSDPIQVKTQQGVPGQPMNLRAEARSETSITLSWSPPRQESIKEYELLFREGDHGREVGRTF
DPTTSYVVEDLKPNTYAFRLAARSPQGLGAFTPVVRQRTLQSIKPKNFVKMIMKTSVLLSWEFPDNYN
SPTPYKIQYNGLTLDVDGRTTKLLITHLKPHTFYNFVLTNRGSSLGGLQQTVAWTA FNLLNGKPSVAPK
PDADGFIMVYLPDGQSPVPVQSYFIVMVPLRKSRRGGQFLPLGSPEDMDLEELIQDISRLQRRSLRHSRQ
LEVPRPYIARFVLPPTFHPGDQKQYGGFDNRGLEPHRYVLFVLAFLQKSEPTFAASPFSDPFQLDNP
DPQPIVDGEEGLIWWIGPVLAVVFIICIVAILLYKNKPKDSKRKDSEPRTKCLLNNADLAPHHPKDPVEM
RRINFQTPGMLSHPIPIADMAEHTERL KANDSLKLSQYYESIDPGQQTWEHSNLEVNKPKNRYANVIA
YDHSRVILQPIEGIMGSDYINANYVDGYRRQNAYIATQGPLPETFGDFWRMVWEQRSATIVMMTRLEEKS
RIKCDQYWPNRGTETYGFIQVTLTDLTIELATFCVRTFSLHKNKGSSEKREVRQFQFTAWPDHGVPEYPTPF
LAFLRRVKT CNPPDAGPIVHCSAGVGRGTCFVIDAMLERIKPEKTVDVYGHVTLMRSQRNYMVQTEDQ
YSFIHEALLEAVGCGNTEVPARSLYAIQKLAQVEPGEHVTGMELEFKRLANSKAHTSRFISANLPCNKF
KNRLVNIMPYESTRVCLQPIRGVEGSDYINASFIDGYRQQKAYIATQGPLAETTEDFWRMLWENNSTIVV
MLTKLREMGREKCHQYWAERSARYQYFVVDPMAEYNMPQYILREFKVTDARDGQSRTVTRQFQFTDWPEQ
GVPKSGEGFIDFIGVHKTKEQFGQDGPISVHCSAGVGRGTVFITLSIVLERMRYEGVVDIFQTVKMLRT
QRPAMVQTEDEYQFCYQAALYLG SFDHYAT

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online >](#)

Tag:	C-Myc/DDK
Predicted MW:	165.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_570923
Locus ID:	5802
UniProt ID:	Q13332 , Q8NHS7 , Q59FX6
RefSeq Size:	6006
Cytogenetics:	19p13.3
RefSeq ORF:	4503
Synonyms:	PTPSIGMA; R-PTP-S; R-PTP-sigma
Summary:	<p>The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP contains an extracellular region, a single transmembrane segment and two tandem intracytoplasmic catalytic domains, and thus represents a receptor-type PTP. The extracellular region of this protein is composed of multiple Ig-like and fibronectin type III-like domains. Studies of the similar gene in mice suggested that this PTP may be involved in cell-cell interaction, primary axonogenesis, and axon guidance during embryogenesis. This PTP has been also implicated in the molecular control of adult nerve repair. Four alternatively spliced transcript variants, which encode distinct proteins, have been reported. [provided by RefSeq, Jul 2008]</p>
Protein Families:	Druggable Genome, Phosphatase, Transmembrane

Product images:

Coomassie blue staining of purified PTPRS protein (Cat# [TP311163]). The protein was produced from HEK293T cells transfected with PTPRS cDNA clone (Cat# [RC211163]) using MegaTran 2.0 (Cat# [TT210002]).