

Product datasheet for **TP311169M**

Interferon alpha 8 (IFNA8) (NM_002170) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human interferon, alpha 8 (IFNA8), 100 µg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC211169 protein sequence
Red=Cloning site **Green**=Tags(s)

MALTFYLLVALVLSYKSFSSLGCDLPQTHSLGNRRALILLAQMRRISPFSCCLKDRHDFEFPQEEFDDKQ
FQKAQAISVLHEMIQQTFNLFSTKDSSAALDETLLEDFYIELDQQLNDLESCVMQEVGVIESPLMYEDSI
LAVRKYFQRITLYLTEKKYSSCAWEVVR AEIMRSFSLINLQKRLKSKE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 21.8 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_002161](#)

Locus ID: 3445

UniProt ID: [P32881](#)


RefSeq Size: 1039



[View online »](#)

Cytogenetics:	9p21.3
RefSeq ORF:	567
Synonyms:	IFN-alphaB
Summary:	Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase. [UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Antigen processing and presentation, Autoimmune thyroid disease, Cytokine-cytokine receptor interaction, Cytosolic DNA-sensing pathway, Jak-STAT signaling pathway, Natural killer cell mediated cytotoxicity, Regulation of autophagy, RIG-I-like receptor signaling pathway, Toll-like receptor signaling pathway

Product images:



Coomassie blue staining of purified IFNA8 protein (Cat# [TP311169]). The protein was produced from HEK293T cells transfected with IFNA8 cDNA clone (Cat# [RC211169]) using MegaTran 2.0 (Cat# [TT210002]).