

Product datasheet for TP311181L

PDE10A (NM_006661) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phosphodiesterase 10A (PDE10A), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211181 representing NM_006661 Red=Cloning site Green=Tags(s)

MRIEERKSQHLTGLTDEKVKAYLSLHPQVLDEFVSESVAETVEKWLKRKNNKSEDESAPKEVSRVYQDTN
MQGVVYELNSYIEQRLDTGGDNQLLLYELSSIIKIATKADGFALYFLGECNNSLCIFTPPGIKEGKPRLI
PAGPITQGTTSAYVAKSRKTLLEDILGDERFPRGTGLESGETRIQSVLCLPIVTAIGDLIGILELYRHW
GKEAFCLSHQEVATANLAWASVAIHQVQVCRGLAKQTELNDFLLDVSKTYFDNIVAIDSLLEHIMYAKN
LVNADRCALFQVDHKNKELYSDLFDIGEEKEGKPVFKKTKAIRFSIEKGIAGQVARTGEVLNIPDAYADP
RFNREVDLYTGYTTRNILCMPIVSRGVSIGVWQMVNKGSAFSDENNFKMFVFCALALHCANMYHR
IRHSECIYRVTMEKLSYHSICTSEEWQGLMQFTLPVRLCKEIELFHFDIGPFENMWPGIFVYMVHRSCGT
SCFELEKLCRFIMSVKKNYRRVPYHNWKHAVTVAHCMYAILQNNHTLFTDLERKGLLIACLCHDLDRGF
SNSYLQKFDHPLAALYSTSTMEQHHFSQTVSILQLEGHNIFSTLSSEYEQVLEIRKAIATDLALYFG
NRKQLEEMYQTGSLNLNQSHRDRVIGLMMTACDLCSVTKLWPVTKLTANDIYAEFWAEGDEMKKLGIQP
IPMMDRDKKDEVPQGLGFYNAVAIPCYTTLTQILPPTPELLKACRDNLSQWEKVRGEETATWISSPSV
AQKAAASED

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	88.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online >](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_006652
Locus ID:	10846
UniProt ID:	Q9Y233 , A0A1B1UZR0
RefSeq Size:	4576
Cytogenetics:	6q27
RefSeq ORF:	2337
Synonyms:	ADSD2; HSPDE10A; IOLOD; LINC00473; PDE10A19
Summary:	The protein encoded by this gene belongs to the cyclic nucleotide phosphodiesterase family. It plays a role in signal transduction by regulating the intracellular concentration of cyclic nucleotides. This protein can hydrolyze both cAMP and cGMP to the corresponding nucleoside 5' monophosphate, but has higher affinity for cAMP, and is more efficient with cAMP as substrate. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Dec 2011]
Protein Families:	Druggable Genome
Protein Pathways:	Progesterone-mediated oocyte maturation, Purine metabolism

Product images:



Coomassie blue staining of purified PDE10A protein (Cat# [TP311181]). The protein was produced from HEK293T cells transfected with PDE10A cDNA clone (Cat# [RC211181]) using MegaTran 2.0 (Cat# [TT210002]).