

Product datasheet for **TP311201M**

HNF1 alpha (HNF1A) (NM_000545) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human HNF1 homeobox A (HNF1A), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211201 representing NM_000545 Red =Cloning site Green =Tags(s)

MVSKLSQLQTELLAALLESGLSKEALLQALGEPGPYLLAGEGPLDKGESCGGGRGELAEPLNGLGETRGS
EDETDDDDGEDFTPPILKELENLSPEEAAHQKAVETLLQEDPWRVAKMVKSYLQQHNIPQREVDTTGLN
QSHLSQHLNKGTPMKTQKRAALYTWYVRKQREVAQQFTHAGQGGLIEEPTGDELPTKKGRRNRFKWGPAS
QQILFQAYERQKNPSKEERETLVEECNRAECIQRGVSPSQAQGLGSNLVTEVRVYNWFANRRKEEAFRHK
LAMDTYSGPPPGPGPALPAHSSPGLPPPALSPOKSVHGVRYGQPATSETAEVPSSSGGPLVTVSTPLHQ
VSPTGLEPSHLLSTEAKLVAAGGPLPPVSTLTALHSLEQTSPGLNQPPQNLIMASLPGVMTIGPGEP
SLGPTFTNTGASTLVIGLASTQAQSVPVINSMGSSLTTLQPVQFSQPLHPSYQQPLMPPVQSHVTQSPFM
ATMAQLQSPHALYSHKPEVAQYHTHTGLLPQTMLITDTTNSALASLTPTKQVFTSDTEASSESGLHTPAS
QATTLHVPSQDPAGIQHLQPAHRLSASPTVSSSLVLYQSSDSSNGQSHLLPSNHSVITFISTQMASSS
Q

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	67.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000536
Locus ID:	6927
UniProt ID:	P20823 , E0YMI7
RefSeq Size:	3249
Cytogenetics:	12q24.31
RefSeq ORF:	1893
Synonyms:	HNF-1A; HNF1; HNF1alpha; HNF4A; IDDM20; LFB1; MODY3; TCF-1; TCF1
Summary:	The protein encoded by this gene is a transcription factor required for the expression of several liver-specific genes. The encoded protein functions as a homodimer and binds to the inverted palindrome 5'-GTTAATNATTAAC-3'. Defects in this gene are a cause of maturity onset diabetes of the young type 3 (MODY3) and also can result in the appearance of hepatic adenomas. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2015]
Protein Families:	Adult stem cells, Druggable Genome, ES Cell Differentiation/IPS, Transcription Factors
Protein Pathways:	Maturity onset diabetes of the young

Product images:



Coomassie blue staining of purified HNF1A protein (Cat# [TP311201]). The protein was produced from HEK293T cells transfected with HNF1A cDNA clone (Cat# [RC211201]) using MegaTran 2.0 (Cat# [TT210002]).