

Product datasheet for TP311245

KIAA1958 (NM_133465) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human KIAA1958 (KIAA1958), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211245 protein sequence Red =Cloning site Green =Tags(s)

MEDCLHTSSENLSKLVSWAHSHGTICSLIPNLKHLLESGSHGNLTAMWGCSAGHAYHWPLTATCRAGSQE
RVCFQDNRSFNSDSPSIIGVPSETQTSPVERYPGRPVKAKLDCNRTRDSCDFSVCSEPELDETVVEEYED
ENTLFDMVCESSVTDESDFEPTQRPQSIARKRPGVWPSSLHSSSQTMVDECSNDVVIKKIKQEIPED
YYIVANAELTGGVDGPALSLTQMAKPKPQTHAGPSCVGSAKLIPHVTSISTELDPHGMSASPSVISRPI
VQKTARVSLASPNRGPPGTHGTNQVAMQMPVSTSHPNKQISIPLSALQLPGQDEQVASEEFLSHLP SQV
SSCEVALSPSVNTEPEVSSSQQPPVAPAITTEATAQCIPAYSTKLNKFPVFNINDDLNDLCTSAVSPNT
TKATRYALNVWRYWCMTNGLKDHTDITKIPAVKLNELLENFYVTVKKS DGSDFLATSLHAIRRLDRILK
NAGVGFSITSSTFSSTKKLKEKLWVLSKAGMSGARSRNIVYFSLSDEEEMWQAGCLGDDSPITLLSTVV
KYNSQYLNMR TLQEHADLMYGDIELLKDPQNQPYFARTDSVKRESRSGSTRVCHGKIYHEHSRGHKQCPY
CLLYKMYIHRPPTQMEAKSPFYLTARKEATDMGSVWYEEQRMGLRSLRGIVPNLAKKVKLENCENFTFV
SFTQVSRRLGSHSCCQ


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	79 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_597722</u>
Locus ID:	158405
UniProt ID:	<u>Q8N8K9</u>
RefSeq Size:	7575
Cytogenetics:	9q32
RefSeq ORF:	2148

Product images:


Coomassie blue staining of purified KIAA1958 protein (Cat# TP311245). The protein was produced from HEK293T cells transfected with KIAA1958 cDNA clone (Cat# [RC211245]) using MegaTran 2.0 (Cat# [TT210002]).