

Product datasheet for TP311245M

KIAA1958 (NM_133465) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human KIAA1958 (KIAA1958), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211245 protein sequence Red =Cloning site Green =Tags(s)

MEDCLHTSSENLSKLVSWAHSHGTICSLIPNLKHLLESGSHGNLTAMWGCSAGHAYHWPLTATCRAGSQ
 E
 RVCFQDNRSFNDSFPSIIGVPSSETQTSPVERYPGRPVKAKLDCNRTRDSCDFSYCSEPELDETVEEYED
 ENTLFDMVCESSVTDESDSFEPQTQRPQSIARKRPGVVPSSLHSSQTQMVDECSNDVIIKKIKQEIPED
 YYIVANAELTGGVDGPALSLTQMAKPKPQTHAGPSCVGSAKLIPHVTSIAISTELDPHGMSASPSVISRPI
 VQKTARVSLASPNRGPPGTHGTNQQVAMQMPVSTSHPNKQISIPLSALQLPGQDEQVASEEFLSHLPSQ
 V
 SSCEVALSPSVNTEPEVSSSQQPPVAPAITTEATAQCIPAYSTKLNKFPVFNINDDLNDLCTSAVSPNT
 TKATRYALNVWRYWCMTNGLKDHTDITKIPAVKLNELLENFYVTVKKSDGSDFLATSLHAIRRGDLRLK
 NAGVGFSITSSTFSSTKKLKEKLWVLSKAGMSGARSNIVYFSLSDEEEMWQAGCLGDDSPITLLSTV
 KYNSQYLNMRTLQEHADLMYGDIELLKDPQNQPYFARTDSVKRESRSGSTRVCHGKIYHEHSRGHKQCP
 Y
 CLLYKMYIHRPPTQMEAKSPFYLTARKEATDMGSVWYEEQRMGLRSLRGIVPNLAKKVKLENCENFTFV
 SFTQVSRRLGSHSCCQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	79 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.


[View online »](#)

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_597722](#)

Locus ID: 158405


UniProt ID: [Q8N8K9](#)

RefSeq Size: 7575

Cytogenetics: 9q32

RefSeq ORF: 2148

Product images:



Coomassie blue staining of purified KIAA1958 protein (Cat# [TP311245]). The protein was produced from HEK293T cells transfected with KIAA1958 cDNA clone (Cat# [RC211245]) using MegaTran 2.0 (Cat# [TT210002]).