

Product datasheet for TP311269L


C20orf79 (SCP2D1) (NM_178483) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 20 open reading frame 79 (C20orf79), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211269 protein sequence Red =Cloning site Green =Tags(s) MWKRSDHQPKIKAEDGPLVGQFEVLGSPVPEPAMPHPLELSEFESFPVFQDIRLHIREVGAQLVKKVNAVF QLDITKNGKTILRWTIDLKNGSGDMYPGPAPRLPADTVFTIPESVFMELVLGKMNPQKAFLAGKFKVSGKV LLSWKLERVFKDWAKF TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	17.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_848578
Locus ID:	140856
UniProt ID:	Q9UJQ7
RefSeq Size:	677


[View online »](#)

Cytogenetics: 20p11.23
RefSeq ORF: 468
Synonyms: C20orf79; dj1068E13.2; HSD22

Product images:


Coomassie blue staining of purified SCP2D1 protein (Cat# [TP311269]). The protein was produced from HEK293T cells transfected with SCP2D1 cDNA clone (Cat# [RC211269]) using MegaTran 2.0 (Cat# [TT210002]).