

Product datasheet for **TP311316L**

PCDHGC5 (NM_032407) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protocadherin gamma subfamily C, 5 (PCDHGC5), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211316 protein sequence Red =Cloning site Green =Tags(s)
	<p>MGPKTLPQLAGKWQVLCMLSLCCWGWVSGQLRYSVVEESEPGLVGNVAQDLGLKMTDLLSRRLQLGSE EE</p> <p>NGRYFSLMSGALAVNQKIDRESLCGASTSCLLPVQVTEHPLELIRVEVEILDNDNSPSFATPEREM RISESAASGARFPLDSAQDPDVGNTNTVSFYTLSPNSHFSLNVKTLKDGKPFPELVLEQQLDREAQARHQL VLTAVDGGTTPARSGTTLISVIVLDINDNAPTFQSSVLRVGIPENAPIGTLRLNATDPDEGTNGQLDYS FGDHTSEAVRNLFGLDPSSGAIHVLGPIDFEESRFYEHARARDQGGPAMEGHCVIQVDVGDVNDNAPEV LLASLANPVLESTPVGTWVGLFNVRDRDSGRNGEVSLDISPDLPFQIKPSENHYSLTSQPLDREATSHY IIELLASDAGSPSLHKHLTIRLNISDVNDNAPRFNQQLYTAYILENRPPGSLCTVAASDPDTGDNARLT YSIVGNQVQGAPASSFVYNPEDGRIFAQRTFDYELLQMLQIVVGVRDSGSPPLHANTSLHVFVLDENDN APAVLHPRPGWEHSAPQRLPRSAPPGSLVTKVTAVDADAGHNAWLSYSLLPQSTAPGLFLVSTHTGEVRT ARALLEDDSDTQQVVLVRDNGDPSLSSTATVLLVLEDEDPEEMPKSSDFLIHPPERSDLTYLIVALAT VSLLSLVTFTFLSAKCLQGNADGDGGGGGCCRRQDPSPDFYKQSSPNLQVSSDGTKYMEVTLRPTDSQ SHCYRTCFSASDGSDFTLRPLSVQQPTALALEPDAIRSRSNLRLERSQVRGSAPPRATPGGGTGEAAR PHKGLNLHPLLSGRLGRWLRSTRFSGSLDRLRGTRVAD</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	92 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_115783
Locus ID:	56097
UniProt ID:	Q9Y5F6
RefSeq Size:	2637
Cytogenetics:	5q31.3
RefSeq ORF:	2634
Synonyms:	PCDH-GAMMA-C5
Summary:	<p>This gene is a member of the protocadherin gamma gene cluster, one of three related clusters tandemly linked on chromosome five. These gene clusters have an immunoglobulin-like organization, suggesting that a novel mechanism may be involved in their regulation and expression. The gamma gene cluster includes 22 genes divided into 3 subfamilies. Subfamily A contains 12 genes, subfamily B contains 7 genes and 2 pseudogenes, and the more distantly related subfamily C contains 3 genes. The tandem array of 22 large, variable region exons are followed by a constant region, containing 3 exons shared by all genes in the cluster. Each variable region exon encodes the extracellular region, which includes 6 cadherin ectodomains and a transmembrane region. The constant region exons encode the common cytoplasmic region. These neural cadherin-like cell adhesion proteins most likely play a critical role in the establishment and function of specific cell-cell connections in the brain. Alternative splicing has been described for the gamma cluster genes. [provided by RefSeq, Jul 2008]</p>
Protein Families:	Transmembrane

Product images:



Coomassie blue staining of purified PCDHGC5 protein (Cat# [TP311316]). The protein was produced from HEK293T cells transfected with PCDHGC5 cDNA clone (Cat# [RC211316]) using MegaTran 2.0 (Cat# [TT210002]).