

## Product datasheet for **TP311346**

### SPRR4 (NM\_173080) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human small proline-rich protein 4 (SPRR4), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211346 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MSSQQQQRQQQQCPPQRAQQQVKQPCPPPVKCQETCAPKTKDPCAPQVKKQCPPKGTIIPAQQKCPSA QQASKSKQK  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	8.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_775103</a>
Locus ID:	163778
UniProt ID:	<a href="#">Q96P11</a>
RefSeq Size:	766
Cytogenetics:	1q21.3




[View online »](#)

RefSeq ORF: 237

Summary: Cross-linked envelope protein of keratinocytes. Involved in UV-induced cornification.  
[UniProtKB/Swiss-Prot Function]

### Product images:



Coomassie blue staining of purified SPRR4 protein (Cat# TP311346). The protein was produced from HEK293T cells transfected with SPRR4 cDNA clone (Cat# [RC211346]) using MegaTran 2.0 (Cat# [TT210002]).