

Product datasheet for TP311389L

DDR2 (NM_006182) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human discoidin domain receptor tyrosine kinase 2 (DDR2), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211389 representing NM_006182 Red=Cloning site Green=Tags(s)

MILIPRMLLVLFLLLPILSSAKAQVNPAICRYPLGMSGGQIPDEDITASSQWSESTAAKYGRLDSEEGDG
AWCPEIPVEPDDLKEFLQIDLHTLHFITLVGTQGRHAGGHGIEFAPMYKINYSRDGTRWISWRNRHGKQV
LDGNSNPYDIFLKDLEPPIVARFVRFPVTDHSMNVCMRVELYGCWLDGLVSYNAPAGQQFVLPGGSII
YLNDSVYDGA VGYSMTEGLGQLTDGVSLDDFTQTHEYHVWPGYDYVGVWRNESATNGYIEIMFEFDRIRN
FTTMKVHCNNMFAKGVKIFKEVQCYFRSEASEWEPNAISFPLVLDDVNPSARFVTVPLHHRMASAIKCQY
HFADTWMMFSEITFQSDAAMYNNSEALPTSPMAPTTYDPMLKVDDSNTRILIGCLVAIIFILLAIIVIL
WRQFWQKMLEKASRRMLDDEMTVLSLPSDSSMFNNRSSSPSEQGSNSTYDRIFLRPDYQEPSRLIRK
LPEFAPGEEESGCSGVVKPVQPSGPEGVPHYAEADIVNLQGVTTGGNTYSVPAVTMDLLSGKDVAVEEFPR
KLLTFKEKLGEGQFGEVHLCVEGMEKFKDKDFALDVSANQPVLVAVKMLRADANKNARNDFLKEIKIMS
RLKDPNIIHLLAVCITDDPLCMITEYMENGDNLNQLSRHEPPNSSSSDVRTVSYTNLKFMATQIASGMKY
LSSLNFVHRDLATRNCLVGKNYTIKIADFGMSRNLYSGDYRIQGRAVLPIRWMSWESILLGKFTTASDV
WAFGVTLWETFTFCQEQPYSQLSDEQVIENTGEFFRDQGRQTYLPQPAICPDSVYKMLMLSCWRRDTKNRP
SFQEIHLLLLLQQGDE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	94.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_006173
Locus ID:	4921
UniProt ID:	Q16832 , A0A024R906
RefSeq Size:	3172
Cytogenetics:	1q23.3
RefSeq ORF:	2565
Synonyms:	MIG20a; NTRKR3; TKT; TYRO10; WRCN
Summary:	This gene encodes a member of the discoidin domain receptor subclass of the receptor tyrosine kinase (RTKs) protein family. RTKs play a key role in the communication of cells with their microenvironment. The encoded protein is a collagen-induced receptor that activates signal transduction pathways involved in cell adhesion, proliferation, and extracellular matrix remodeling. This protein is expressed in numerous cell types and may also be involved in wound repair and regulate tumor growth and invasiveness. Mutations in this gene are the cause of short limb-hand type spondylometaepiphyseal dysplasia. [provided by RefSeq, Aug 2017]
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane

Product images:



Coomassie blue staining of purified DDR2 protein (Cat# [TP311389]). The protein was produced from HEK293T cells transfected with DDR2 cDNA clone (Cat# [RC211389]) using MegaTran 2.0 (Cat# [TT210002]).