

## Product datasheet for TP311516L


### C8orf74 (NM\_001040032) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 8 open reading frame 74 (C8orf74), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211516 protein sequence <span style="color: red;">Red</span> =Cloning site <span style="color: green;">Green</span> =Tags(s)  MALLTPQGVKEVFQFQRPQGRERLRLLNWEEFDEQRDSRRSILLDTLYESIIFAVGKGFVVEVAQVVK FTEELLRETKGCSITEAVTILGNKLRDYRGHFNTTHLLALCDYFHHTFIRHYKLYQVVLGQDQQVDLTVA HLEVCMPHPLPLAEGMDRDLWIHEQQVATLTEAEAKRADVLLLKEALRLERENSLQKAFAAAAAPAQP G QVLERQELESICQAVHTQMELLQELLQRQIQNTFAILDCLKKTLNLNAPTPIPPITSHAGQEEALK PQRASKGKKAKARK  <span style="color: red;">TR</span> <span style="color: green;">TRPLEQKLISEEDLA</span> <span style="color: green;">ANDILDYKDDDDKV</span>
Tag:	C-Myc/DDK
Predicted MW:	33.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001035121</a></u>
Locus ID:	203076


[View online »](#)

UniProt ID: Q6P047  
 RefSeq Size: 1080  
 Cytogenetics: 8p23.1  
 RefSeq ORF: 882

**Product images:**


Coomassie blue staining of purified C8orf74 protein (Cat# [TP311516]). The protein was produced from HEK293T cells transfected with C8orf74 cDNA clone (Cat# [RC211516]) using MegaTran 2.0 (Cat# [TT210002]).