

Product datasheet for TP311711M

SLURP2 (NM_177458) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human Ly6/neurotoxin 1 (LYNX1), transcript variant SLURP2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211711 representing NM_177458 <div> <div>Red</div>=Cloning site <div>Green</div>=Tags(s) </div> <p> MQLGTGLLLAAVLSLQLAAAEAIWCHQCTGFGGCSHGSRCLRDSTHCVTATRVLSNTEDLPLVTKMCHI GCPDIPSLGLGPYVSIACCQTSLCNHD <div> <div>TR</div> <div>TRPLEQKLISEEDLAANDILDYKDDDDKV</div> </div> </p>
Tag:	C-Myc/DDK
Predicted MW:	10 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_803253
Locus ID:	432355
UniProt ID:	P0DP57
RefSeq Size:	589
Cytogenetics:	8q24.3



[View online »](#)

RefSeq ORF: 291

Synonyms: SLURP-2

Summary: This gene encodes a novel, secreted member of the Ly6/uPAR (LU) superfamily of proteins containing the unique three-finger LU domain. This gene is mainly expressed in epithelial cells, including skin and keratinocytes, and is up-regulated in psoriatic skin lesions, suggesting its involvement in the pathophysiology of psoriasis. Alternatively spliced transcript variants have been found for this gene. Read-through transcription from the neighboring upstream gene (LYNX1) generates naturally-occurring transcripts (LYNX1-SLURP2) that encode a fusion protein comprised of sequence sharing identity with each individual gene product. [provided by RefSeq, Sep 2017]

Product images:



Coomassie blue staining of purified SLURP2 protein (Cat# [TP311711]). The protein was produced from HEK293T cells transfected with SLURP2 cDNA clone (Cat# [RC211711]) using MegaTran 2.0 (Cat# [TT210002]).