

Product datasheet for **TP311725L**

KIR2DS2 (NM_012312) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant protein of human killer cell immunoglobulin-like receptor, two domains, short cytoplasmic tail, 2 (KIR2DS2), 1 mg |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >Peptide sequence encoded by RC211725 Blue =ORF Red =Cloning site Green =Tag(s) |

MSLTVSMACVGFFLLQGAWPHEGVHRKPSLLAHPGPLVKSEETVILQCWSDVRFELHLLHREGKYKDT
LHLIGEHHHDGVSKANFSIGPMMQDLAGTYRCYGSVTHSPYQLSAPSDPLDIVITGLYEKPSLSAQPGPT
VLAGESVTLSCSRSSYDMYHLSREGEAHERRFSAGPKVNGTFQADFPLGPATHGGTYRCFGSFRDSPY
EWSNSSDPLLVSVTGNPSNSWPSPTPESSKTGNPRHLHVLIGTSVKIPFTILLFLLHRWCNKKNA
VMDQEPAGNRTVNSEDSDEQDHQEVSYA
TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Recombinant protein using RC211725 also available, [TP311725M](#)


| | |
|----------------|--|
| Tag: | C-Myc/DDK |
| Predicted MW: | 31.5 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Preparation: | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps. |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_036444 |



[View online »](#)

| | |
|---------------|--|
| Locus ID: | 100132285 |
| UniProt ID: | P43631 |
| RefSeq Size: | 1573 |
| Cytogenetics: | 19q13.4 |
| RefSeq ORF: | 912 |
| Synonyms: | 183Act1; CD158b; CD158J; cl-49; KIR-2DS2; KIR2DL1; NKAT-5; NKAT5 |
| Summary: | <p>Killer cell immunoglobulin-like receptors (KIRs) are transmembrane glycoproteins expressed by natural killer cells and subsets of T cells. The KIR genes are polymorphic and highly homologous and they are found in a cluster on chromosome 19q13.4 within the 1 Mb leukocyte receptor complex (LRC). The gene content of the KIR gene cluster varies among haplotypes, although several "framework" genes are found in all haplotypes (KIR3DL3, KIR3DP1, KIR3DL4, KIR3DL2). The KIR proteins are classified by the number of extracellular immunoglobulin domains (2D or 3D) and by whether they have a long (L) or short (S) cytoplasmic domain. KIR proteins with the long cytoplasmic domain transduce inhibitory signals upon ligand binding via an immune tyrosine-based inhibitory motif (ITIM), while KIR proteins with the short cytoplasmic domain lack the ITIM motif and instead associate with the TYRO protein tyrosine kinase binding protein to transduce activating signals. The ligands for several KIR proteins are subsets of HLA class I molecules; thus, KIR proteins are thought to play an important role in regulation of the immune response. This gene represents a haplotype-specific family member that encodes a protein with a short cytoplasmic tail. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]</p> |

Product images:



Coomassie blue staining of purified KIR2DS2 protein (Cat# [TP311725]). The protein was produced from HEK293T cells transfected with KIR2DS2 cDNA clone (Cat# [RC211725]) using MegaTran 2.0 (Cat# [TT210002]).