

Product datasheet for **TP311728L**

CYP11A1 (NM_001099773) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens cytochrome P450, family 11, subfamily A, polypeptide 1 (CYP11A1), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211728 representing NM_001099773 Red =Cloning site Green =Tags(s) MAPEATKNFLPLLDVSRDFVSVLHRRRIKAGSGNYSGLDISDDLFRFAFESITNVIFGERQGMLEEVNP EAQRFIDAIYQMFHTSVPMNLNPPDLFRLFRTKTWKDHVAAWDVIFSKADIYTQNFYWELRQKGSVHHDY RGILYRLLGDSKMSFEDIKANVTEMLAGGVDTTSMTLQWHLYEMARNLKVQDMLRAEVLAARHQAQGDMA TMLQLVPLLKASIKETLRLHPISVTLQRYLVNDLVLRDYMIPAKTLVQVAIYALGREPTFFFDPENFDPT RWLSKDKNITYFRNLGFGWGVQRQCLGRRIAELEMTIFLINMLENFRVEIQHLSDVGTTFNLILMPEKPIS FTFWPFNQEATQQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	42 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001093243</u>



[View online »](#)

Locus ID: 1583

UniProt ID: [P05108](#), [P05108-2](#)

RefSeq Size: 2010

Cytogenetics: 15q24.1

RefSeq ORF: 1089


Synonyms: CYP11A1; CYPXIA1; P450SCC

Summary: This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the mitochondrial inner membrane and catalyzes the conversion of cholesterol to pregnenolone, the first and rate-limiting step in the synthesis of the steroid hormones. Two transcript variants encoding different isoforms have been found for this gene. The cellular location of the smaller isoform is unclear since it lacks the mitochondrial-targeting transit peptide. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, P450

Protein Pathways: C21-Steroid hormone metabolism, Metabolic pathways

Product images:



Coomassie blue staining of purified CYP11A1 protein (Cat# [TP311728]). The protein was produced from HEK293T cells transfected with CYP11A1 cDNA clone (Cat# [RC211728]) using MegaTran 2.0 (Cat# [TT210002]).