

Product datasheet for **TP311745L**

SIN1 (MAPKAP1) (NM_001006617) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human mitogen-activated protein kinase associated protein 1 (MAPKAP1), transcript variant 1, full length, with C-terminal MYC/DDK tag, expressed in HEK293 cells, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211745 representing NM_001006617 Red =Cloning site Green =Tags(s)

MAFLDNPTIILAHIRQSHVTSDDTGMCEMVLIDHDVDLEKIHPPSMPGDSGSEIQGSNGETQGYVYAQSV
DITSSWDFGIRRRSNTAQRLELRKERQNIQCKNIQWKERNKQSAQELKSLFEKKSLKEKPPISGKQS
ILSVRLEQCPLQLNPNFNEYSKFDGKGVGTATTAKKIDVYLPVHSSQDRLLPMTVVTMASARVQDLIGLI
CWQYTSEGREPKLNDNVSAAYCLHIAEDDGEVDTDFPPLDSNEPIHKFGFSTLALVEKYSSPGLTSKESLF
VRINAAHGFSLIQVDNNTKVTMKEILLKAVKRRKGSQKVGSPQYRLEKQSEPNVAVDLSTLESQSAWEFC
LVRENSSRADGVFEEDSQIDIATVQDMLSSHYYKSFVSMIHRLRFTTQVQLGISGDKVEIDPVTNQKAS
TKFWIKQKPISIDSDLLCACDLAEKSPSHAIFKLTLYLSNHDKHLYFESDAATVNEIVLKVNYILESRA
STARADYFAQKQRKLNRRTSFSFQKEKKSQQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	58.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_001006618](#)

Locus ID: 79109

UniProt ID: [Q9BPZ7](#)

RefSeq Size: 3395

Cytogenetics: 9q33.3

RefSeq ORF: 1566

Synonyms: JC310; MIP1; SIN1; SIN1b; SIN1g

Summary: This gene encodes a protein that is highly similar to the yeast SIN1 protein, a stress-activated protein kinase. Alternatively spliced transcript variants encoding distinct isoforms have been described. Alternate polyadenylation sites as well as alternate 3' UTRs have been identified for transcripts of this gene. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified MAPKAP1 protein (Cat# [TP311745]). The protein was produced from HEK293T cells transfected with MAPKAP1 cDNA clone (Cat# [RC211745]) using MegaTran 2.0 (Cat# [TT210002]).