

Product datasheet for **TP311834L**

SPATA5 (NM_145207) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human spermatogenesis associated 5 (SPATA5), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211834 representing NM_145207 Red =Cloning site Green =Tags(s)

MSSKKNRRLNQSAENGSSLSAASSCAEARAPSAGSDFAAATSGTLTVTNLLEKVDDKIPKTFQNSLIHL
GLNTMKSANICIGRPVLLTSLNGKQEVYTAWPMAGFPGGKVGLSEMAQKNVGVPRGDIAIQVQPLVGAVLQ
AEEMDVALSDKDMEINEEELTGCILRKLKGKIVLPGNFLYCTFYGRPYKLQVLRVKGADGMILGGPQSDS
DTDAQRMFAEQSSMETSSLESLQLSLDLEDTQIPTSRSTPYKPIDDRITNKASDVLLDVTQSPGDGSG
LMLEEVTGLKCNFESAREGNEQLTEERLLKFSIGAKCNTDTFYFISSTTRVNFTEIDKNSKEQDNQFKV
TYDMIGGLSSQLKAIREIIEPLKQPELFSYGIAPRGVLLYGPPGTGKTMARAVANEVGAIVSVING
PEIISKFYGETEAKLRQIFAEATLRHPSIIFIDELDALCPKREGAQNEVEKRVASLLTMDGIGSEVSE
GQVLVLGATNRPHALDAALRRPGRFDKEIEIGVNAQDRDLILQKLLRRVPHLLTEALLQLANSAHGYY
GADLKVLCNEAGLCALRRILKKQPNLPDVKVAGLVKITLKDFLQAMNDIRPSAMREIAIDVPNVSWSDIG
GLESIKLLEQAVEWPLKHPESFIRMGIQPPKGVLLYGPPGCSKTMIKALANESGLNFLAIKPELMNK
YVGESERAVRETFRKARAVAPSIIFDELDAVERGSSLGAGNVADRVLAQLLEMDGIEQLKDVTILA
ATNRPDRIDKALMRPGRIDRIIYVPLPDAATRREIFKLQFHSMVPSNEVDLDELILQTDAYS GA EIVAVC
REAALLALEEDIQANLIMKRHFTQALSTVTPRIPESLRRFYEDYQEKSGLHTL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	97.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_660208
Locus ID:	166378
UniProt ID:	Q8NB90
RefSeq Size:	3359
Cytogenetics:	4q28.1
RefSeq ORF:	2679
Synonyms:	AFG2; EHLMRS; SPAF
Summary:	This gene encodes a member of the ATPase associated with diverse activities family, whose members are defined by a highly conserved ATPase domain. Members of this family participate in diverse cellular processes that include membrane fusion, DNA replication, microtubule severing, and protein degradation. The protein encoded by this gene has a putative mitochondrial targeting sequence and has been proposed to function in maintenance of mitochondrial function and integrity during mouse spermatogenesis. Allelic variants in this gene have been associated with epilepsy, hearing loss, and cognitive disability syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016]

Product images:



Coomassie blue staining of purified SPATA5 protein (Cat# [TP311834]). The protein was produced from HEK293T cells transfected with SPATA5 cDNA clone (Cat# [RC211834]) using MegaTran 2.0 (Cat# [TT210002]).