

## Product datasheet for TP311843L

### CCDC43 (NM\_001099225) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 43 (CCDC43), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211843 representing NM_001099225 Red=Cloning site Green=Tags(s)

MAAPSEVAAIAPGEGDGGGGFGSWLDGRLEALGVDRVYVYAYILGILQEEEEEEKLDALQGILSAFLEE  
DSLLNICKEIVERWSETQNVVTKVKKEDDEVQAIATLIEKQAQIVVKPRMVSEEEKQRKAALLAQYADVTD  
EEDSVPKHQCGRCP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	16.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001092695</a>
Locus ID:	124808
UniProt ID:	<a href="#">Q96MW1</a>




[View online »](#)

RefSeq Size: 2088

Cytogenetics: 17q21.31

RefSeq ORF: 462

### Product images:



Coomassie blue staining of purified CCDC43 protein (Cat# [TP311843]). The protein was produced from HEK293T cells transfected with CCDC43 cDNA clone (Cat# [RC211843]) using MegaTran 2.0 (Cat# [TT210002]).