

Product datasheet for **TP311875L**

Clusterin (CLU) (NM_001831) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human clusterin (CLU), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC211875 representing NM_001831 Red =Cloning site Green =Tags(s)

MQVCSQPQRGCVREQSAINTAPPSAHNAASPGGARGHRVPLTEACKDSRIGGMMKTLFFVGLLLTWESG
QVLGDQTVSDNELQEMSNQGSKYVNKEIQNAVNGVKQIKTLIEKTNEERKTLLSNLEEAKKKKEDALNET
RESETKLKELPGVCNETMMALWEECKPCLKQTCMKFYARVCRSGSLVGRQLEEFNLQSSPFYFWMNGDR
IDSLLENDRQQTHMLDVMQDHFSSRASSIIDELFQDRFFTRPQDTYHYLPFSLPHRRPHFFFPKSRIVRS
LMPFSPYEPLNFHAMFQPFLEMIHEAQQAMDIHFHSPAFQHPPTFEFIREGDDDRTVCREIRHNSTGCLRM
KDQCDKCREILSVCSTNNPSQAKLRRELDLQVAERLTRKYNELLSYQWKMLNTSSLLEQLNEQFNW
VSRLANLTQGEDQYYLRVTTVASHTSDSDVPSGVTEVVVKLFSDPITVTPVEVSRKNPKFMETVAEKA
LQEYRKKHREE

SGPTRRRLE**QKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	57.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	Probe target (PMID: 28103719)
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_001822](#)

Locus ID: 1191

UniProt ID: [P10909](#)

RefSeq Size: 2859

Cytogenetics: 8p21.1


RefSeq ORF: 1503

Synonyms: AAG4; APO-J; APOJ; CLI; CLU1; CLU2; KUB1; NA1/NA2; SGP-2; SGP2; SP-40; TRPM-2; TRPM2


Summary: The protein encoded by this gene is a secreted chaperone that can under some stress conditions also be found in the cell cytosol. It has been suggested to be involved in several basic biological events such as cell death, tumor progression, and neurodegenerative disorders. Alternate splicing results in both coding and non-coding variants.[provided by RefSeq, May 2011]

Protein Families: Druggable Genome, Secreted Protein

Product images:



Schematic representation of the conjugation reaction between PEG-functionalized liposomes and the protein CLU via a carboxy-to-amine reaction in the presence of the crosslinker and stabilizer EDC and NHS. CLU Protein CLU from OriGene ([TP311875]) was used in the conjugation reaction. Figure cited from J Liposome Res, PMID: 28103719



Coomassie blue staining of purified CLU protein (Cat# [TP311875]). The protein was produced from HEK293T cells transfected with CLU cDNA clone (Cat# [RC211875]) using MegaTran 2.0 (Cat# [TT210002]).